





# Nanopore sequencing provides superior MGMT promoter methylation evaluation to conventional techniques

Skarphedinn Halldorsson  <sup>\*</sup>1, Richard Nagymihaly<sup>1</sup>, Areeba Patel <sup>2</sup>, Petter Brandal<sup>3</sup>,  
Ioannis Panagopoulos<sup>3</sup>, Felix Sahm <sup>2</sup>, and Einar Vik-Mo  <sup>†</sup>1

<sup>1</sup>Vilhelm Magnus Laboratory, Institute for Surgical Research, Oslo University Hospital

<sup>2</sup>Department of Neuropathology, University Hospital Heidelberg, Heidelberg, Germany

<sup>3</sup>Section for Cancer Cytogenetics, Institute for Cancer Genetics and Informatics, Oslo University Hospital

## Abstract

**Rationale:** Resistance of glioblastoma to the alkylating agent temozolomide may result from the expression of the DNA repair enzyme O6-methylguanine-DNA methyltransferase (MGMT). Methylation of the MGMT promoter region has been correlated with responsiveness to temozolomide, but there is no consensus on the most accurate method to determine methylation. Conventional methods have limitations such as the need for bisulphate treatment and amplification. Long-read Nanopore sequencing offers methylation analysis of native DNA without the need for bisulphate treatment or amplification. Combined with recent advancements in targeting methods, it provides a modern, cost-effective approach to MGMT promoter methylation analysis.

**Methods:** In this study, we analyzed 148 CNS tumors using Nanopore sequencing and compared the results to data obtained using pyrosequencing or methylation bead arrays. We used ONT MinION flow cells to run single or barcoded (multiplex) assays, following a CRISPR/Cas9 protocol, and included results from adaptive sequencing runs. We then compared the methylation data to results from existing methods.

**Results:** We found a 92% correlation between pyrosequencing of 4 CpGs in the MGMT promoter and nanopore sequencing. We observed that samples were liable to be overestimated with ONT, especially with low methylation (<10%) status. We could also re-create classification by the MGMT STP27 algorithm with data from nanopore sequencing. However, using ONT, we were able to include an additional 94 CpGs of the MGMT promoter

---

\*skarphedinn.halldorsson@rr-research.no

†einar.vik-mo@rr-research.no

region in the analysis. Data clustering revealed a robust difference between unmethylated and methylated samples that could be used for patient stratification.

**Discussion:** Our findings demonstrate that ONT is a capable method for replacing pyrosequencing, or methylation bead-array, providing high-confidence results within a few hours of sequencing. The extension of the analysis to the 98 CpG islands in the MGMT promoter region enables further exploration of the correlation between methylation status and additional clinical parameters. However, improving ONT protocols and methodology is necessary to fully replace pyrosequencing in a routine setting.

**Keywords:** MGMT promoter methylation, Nanopore sequencing, CRISPR/Cas9, Glioblastoma

## Introduction

Glioblastoma multiforme (GBM) is a highly aggressive and deadly form of brain cancer, characterized by rapid growth, invasiveness and resistance to treatments. It is the most common and most aggressive type of primary malignant brain tumor in adults [Ostrom et al., 2020]. Despite advances in treatment options, the prognosis for newly diagnosed GBM patients remains poor, with a median survival of less than 15 months [Stupp et al., 2017]. Standard treatment for GBM involves surgical resection of the tumor followed by a combination radiation and chemotherapy. A promising approach to treating GBM involves the use of temozolomide (TMZ), a chemotherapy drug that has been shown to extend the lives of patients when used in combination with radiation therapy [Stupp et al., 2009]. TMZ is an alkylating agent that induces DNA damage by methylation of O-6 guanine residues in dividing cells, leading to DNA damage and apoptosis [Zhang et al., 2011]. However, the effects of TMZ are countered by the DNA repair enzyme O-6-methylguanine DNA methyltransferase (MGMT). MGMT expression is regulated via methylation of the promoter region. The presence of MGMT promoter methylation has been associated with increased survival in glioblastoma patients treated with temozolomide and radiation therapy [Hegi et al., 2019]. Methylation of the MGMT promoter is believed to silence expression, thereby increasing sensitivity of GBM tumor cells to TMZ. MGMT promoter methylation is therefore an important prognostic factor for the management and treatment of GBM [Christmann et al., 2011]. Despite its potential benefits, TMZ can cause a range of side effects, including nausea, vomiting, fatigue, and low blood-cell counts. More severe side-effects such as blood-clots, seizures and liver damage have also been reported. TMZ administration should

be limited to patients that may truly benefit from it and withheld from patients that most likely will only experience the side effects without any improvement in survival.

There is currently no consensus on the optimal method to determine MGMT promoter methylation. Testing is typically performed on tumor tissue samples using techniques such as methylation-specific PCR (MSP), pyrosequencing (PSQ) or methylation bead array. These methods all rely on bisulphate conversion of native tumor DNA prior to analysis and only include a fraction of the 98 CpG sites reported in the promoter/enhancer region of MGMT [Johannessen et al., 2018].

In recent years, advances in sequencing technology have allowed for more sensitive and accurate detection of DNA methylation. Nanopore sequencing, which uses a nanopore-based sensor to detect changes in electrical current as DNA molecules pass through the pore has the ability to detect modifications to the DNA molecule, such as methylation, directly from the signal [Jain et al., 2016]. Due to the long-read nature of nanopore sequencing, it also affords methylation analysis of far longer sequences than either MSP or pyrosequencing. Consequently, nanopore sequencing offers an overview of all the CpG islands of the MGMT promoter region, using native genomic DNA without manipulation which can be both time and cost efficient in a clinical setting [Laver et al., 2015].

Here, we present the results of nanopore sequencing of the promoter region of the MGMT gene in 148 CNS tumors, including 91 GBMs. Results were produced either by CRISPR/Cas9 targeted sequencing of the MGMT promoter region [Wongsurawat et al., 2020] or as part of an adaptive sampling panel [Patel et al., 2022]. We show that nanopore sequencing of the MGMT promoter region, even at low sequencing depth, can accurately recreate the results of pyrosequencing or Illumina 450K bead array. Unsupervised clustering of samples based on methylation of all 98 CpG sites in the MGMT promoter indicates the presence of subgroups within both methylated and unmethylated samples of unknown clinical significance.

## Materials and Methods

### Patients and samples

Samples from three independent cohorts were included into this study; 1) Retrospective analysis of DNA from 68 CNS tumor entities provided by the Institute for Cancer Genetics and Informatics, Oslo University Hospital that had previously been analysed for MGMT promoter methylation

with the Qiagen<sup>®</sup> MGMT pyrosequencing kit. 2) Retrospective analysis of 67 sequences generated as part of the *Rapid-CNS* adaptive sampling pipeline [Patel et al., 2022] that had previously been analysed by Illumina<sup>®</sup> methylation 450K bead array. 3) DNA extracted from 16 primary glioma biopsies that were operated at Oslo University Hospital. A separate biopsy was analysed with the Qiagen<sup>®</sup> MGMT pyrosequencing kit at the [Neuropathology SOMETHING]. Table 1 provides an overview of samples used in this study.

## Sample preparation and Nanopore sequencing

Between 10 and 25 mg of tissue were used to extract genomic DNA (Merck's GenElute<sup>™</sup> Mammalian Genomic DNA Miniprep kit) following the manufacturer's protocol. Purity and concentration of DNA samples was determined using NanoDrop<sup>™</sup> One and Qubit<sup>™</sup> 4 Fluorometers (Thermo Fischer Scientific). Isolated DNA was stored at -20°C until analysis. Cas9 mediated targeted sequencing was performed with the Cas9 Sequencing Kit (Oxford Nanopore Technologies) according to the manufacturers protocol (version ENR\_9084\_v109\_revR\_04Dec2018). Briefly, Cas9 ribonucleoprotein complexes (RNPs) were created by mixing equimolar concentrations (100 µM) of crRNA and trans-activating elements (tracrRNA) to HiFi<sup>®</sup> Cas9 enzyme (IDT). Dephosphorylated gDNA (2-5 µg) was cleaved and dA-tailed with Cas9 RNPs and Taq polymerase. Finally, sequencing adaptors were ligated to cleaved fragments and the final DNA library cleaned with AMPure XP beads (Beckman Coulter). Barcodes were applied to a number of samples to allow multiplexing of five samples based on an experimental protocol from Oxford Nanopore Technologies. Purified DNA libraries were loaded onto R9.4.1. flow cells on MinION Mk1B or Mk1C devices and sequenced for 4-24 hours. Individual flowcells were flushed and re-used up to four times for single samples and twice for multiplexed samples. A minimum pore-count of 300 was deemed sufficient for a single sample, 800 for multiplexed samples. Raw fast5 sequences of all fragments mapping to the MGMT promoter in the Rapid-CNS data were provided for re-analysis.

## Primers

All primers were purchased from Integrated DNA Technologies, IDT (Leuven, Belgium). Previously published primers were used to target the MGMT promoter [Wongsurawat et al., 2020], termed MGMT-left-1 (GCCAACCACGTTAGAGACAATGG); MGMT-left-2 (ATGAGGGGCCCCACTAATTGA); MGMT-right-1 (CCGTAATGTCGGTTATAACACCG) and MGMT-right-2 (GTACGGAGCTATACTCAGGT),

which yielded a fragment of 2,522 bp. In order to reduce background noise from sequencing and expand the size of the fragment, a new pair of reverse primers was designed, MGMT-right-3 (CTGGAATCGCATTCCAGTAGTGG) and MGMT-right-4 (ACTTCGCAAGCATCACAGGTAGG) providing a fragment of 4,800 bp.

## Data analysis

Raw sequences were basecalled, methylation called and mapped (hg19, chromosome 10) using the Megalodon toolbox (version 2.5.0 built on guppy version 6.2.7) from Oxford Nanopore Technologies (<https://github.com/nanoporetech/megalodon>). Methylation percentages of individual CpG sites were compiled using custom scripts in R. All statistical analysis was performed in R.

## Results

A total of 153 samples from 148 patients were analyzed for MGMT promoter methylation, consisting of 91 GBM samples, 23 IDH-glioma samples, and 12 meningioma samples (Figure 2a). Two methods were used to enrich for the region of interest: Cas9 targeted sequencing and adaptive sampling. Cas9 targeted sequencing was applied to 86 samples, 46 of which were run as single samples and 40 that were run as multiplexed groups of 5. 67 samples were analyzed as part of an adaptive sampling pipeline.

Sequence depth of the MGMT promoter region in the samples varied based on method, sequencing time, and DNA and flow-cell quality. Single sample runs produced on average more sequences (mean = 92.1, median = 33) than barcoded runs (mean = 17.2, median = 12) and adaptive sampling (mean = 18.7, median = 15) (Figure 2b). No bias in sequencing depth was observed between methylated and unmethylated samples across Cas9 targeted samples, either single or multiplexed. However, a slight but statistically significant difference in sequence depth was observed between methylated (mean = 23.6) and unmethylated (mean = 16.5) samples created by adaptive sampling ( $p=0.021$ ).

## Nanopore Sequencing versus Pyrosequencing of the MGMT Promoter

Pyrosequencing is a commonly used method to detect MGMT promoter methylation in clinical samples. The Qiagen® MGMT pyrosequencing kit, which detects methylation on 4 CpG sites

(76-79) on the MGMT promoter CpG island, is a common choice in the clinical setting. However, there is no clear consensus on the best cut-off point to classify clinically relevant methylated or unmethylated samples [Brandner et al., 2021]. Cut-off values from 7% to 30% average methylation on the 4 CpGs have been reported (Table 2). Oslo University Hospital considers an average methylation above 10% to be methylated. A subset of our samples (n=68) were initially analyzed using the Qiagen® MGMT pyrosequencing kit before undergoing nanopore sequencing. This allowed us to directly compare the results of the MGMT pyro kit with those of the nanopore sequencing covering the same 4 CpG sites (Figure 3). The correlation between the methylation values of each overlapping CpG site between nanopore and pyrosequencing ranged from 0.78 to 0.88 (Figure 3a). However, the correlation increased to 0.92 when methylation values were averaged across the four CpG sites (Figure 3b).

A 10% average methylation threshold of CpGs 76-79 was applied to the nanopore data to classify MGMT methylated versus unmethylated samples. When these classification results were compared to the classification obtained via pyrosequencing (Figure 3c, left), we found a 91% concordance rate between the two methods (62 out of 68 samples) (Figure 3d, upper). Notably, samples where nanopore sequencing and pyrosequencing gave discordant results were in all cases false positive, classified as methylated by nanopore sequencing but unmethylated by pyrosequencing .

When we applied the same 10% methylation threshold to samples from the adaptive sequencing panel that were previously classified by Illumina® methylation 450K bead array, the concordance between classification methods dropped to 86% (Figure 3c, right). Discordant cases between nanopore sequencing and bead array were both false positives and false negatives (Figure 3d, lower). It's worth noting that these results only apply to the 4 CpGs sequenced by the Qiagen® MGMT pyrosequencing kit.

Illumina® Human Methylation BeadChips (HM-27K, HM-450K, and HM-850K) are microarray-based platforms used to investigate DNA methylation patterns in human tumor samples. Despite detecting the methylation status of tens to hundreds of thousands of CpG sites, these platforms only cover a fraction of the approximately 30 million CpG sites in the human genome. To predict the clinically relevant methylation status of the MGMT promoter, a regression model called *MGMT STP-27* has been developed. This model uses the methylation status of two CpG sites, cg12434587 and cg12981137, as reported by [Bady et al., 2012, Bady et al., 2016].

In the Rapid-CNS study, samples were analyzed by methylation bead array before nanopore sequencing, and the ground truth for MGMT promoter methylation status was inferred from EPIC array results. Methylation values for the two CpG sites represented in the MGMT-STP27 algorithm were extracted from the nanopore data and plotted against each other (Figure 4a). The samples from the Rapid-CNS cohort showed a clear separation between methylated and unmethylated samples based on the methylation percentages of cg12434587 and cg12981137 (Figure 4a, right). In contrast, the samples from the Radium cohort, which were classified as methylated or unmethylated by pyrosequencing, did not show as clear a distinction regarding methylation of the STP27 sites (Figure 4a, left).

Siller *et. al* recently proposed a method for GBM patient stratification by counting the methylation of the 25 CpG sites of the second differentially methylated region (DMR2) in the MGMT promoter using Sanger bisulfite sequencing. Their results showed that the number of methylated CpG sites in DMR2 correlated with patient survival [Siller et al., 2021]. Since nanopore sequencing of the MGMT promoter provides methylation percentages for every individual CpG site within the promoter region and beyond, it is possible to replicate such methods. The results of nanopore sequencing were binarized by applying a methylation cut-off of 10% to each CpG site ( $\geq 10\%$  methylated = methylated,  $< 10\%$  methylated = unmethylated) and summarizing the counts in DMR2. Figure 4b shows a nearly complete separation of methylated and unmethylated samples at  $\geq 15$  methylated CpG sites.

In summary, nanopore sequencing can reproduce the classification results of pyrosequencing, methylation bead array, and bisulfite Sanger sequencing. However, the accuracy of classification depends on the method and cut-offs used to generate the ground truth.

## Unsupervised clustering of samples based on nanopore sequencing

Although classification by bisulphite sequencing methods can be recreated to a reasonable degree with nanopore sequencing data, this does not take advantage of other CpG sites within the designated MGMT promoter CpG island or its shelves and shores that may prove to be relevant for MGMT gene expression. To investigate the impact of methylation at CpG sites not covered by previous methods, we performed hierarchical clustering of 98 CpG sites on the CpG island and included 7 CpGs upstream and 11 CpGs downstream of the CpG island. Unsupervised hierarchical clustering using Ward's method reveals two main clusters that largely correspond to

the classification into methylated and unmethylated samples by pyrosequencing or methylation bead array (Figure 5).

Unmethylated samples exhibit low methylation levels throughout the CpG island, except for the first 5 CpG sites, which are often methylated. On the other hand, methylated samples show a larger gradient of methylation, with higher levels towards either end of the CpG island. This is further supported by the average methylation percentage of each CpG site in methylated and unmethylated samples (Figure 6), which reveals the biggest differences in methylation occur in CpGs 6 through 15 and 71 through 90.

While five samples previously classified as methylated cluster with the otherwise unmethylated samples, one unmethylated sample clusters with methylated samples. This pattern of separation is also evident when unsupervised clustering is performed on GBM samples only (Figure 7). In addition to the robust separation of samples into clusters that largely correspond to the pre-determined methylation status, k-means clustering showed separation of samples in the methylated cluster (Figure 8). Of the 22 samples that cluster with methylated samples, 9 samples fall within what can be described as "very high methylation" cluster. The functional significance of these clusters remains to be determined.

## Survival Analysis

The methylation status of the MGMT promoter is a well-known predictive factor for the overall and progression-free survival of GBM patients receiving Temozolamide treatment [Dovek et al., 2019]. While nanopore methylation profiles were often in agreement with bisulphite sequencing methods, discrepancies were also observed (Figure 5). Therefore, we investigated whether clustering by nanopore sequencing was as effective as the MGMT-pyro kit or EPIC-array for survival prediction. We conducted cas9-targeted nanopore sequencing on 16 additional samples that were simultaneously analyzed by pyrosequencing. In total, we performed survival analysis on 25 primary IDH negative GBM patients (11 females, average age 58.4 years and 14 males, average age 62.7 years) where biopsies were classified by both MGMT-pyro kit and cas9-targeted nanopore sequencing (Table 3). As expected, Kaplan-Meier survival analysis of patients based on pyrosequencing showed a significantly longer overall survival in patients classified as "Methylated" (Figure 9a,  $p=0.0078$ ). Notably, when patients were classified according to unsupervised clustering by nanopore sequencing (Figure 9b), significantly longer survival was observed in "cluster 2" patients ( $p=0.039$ ). Although the sample size is small, our results suggest that classify-



239 ing patients via nanopore sequencing is equally reliable as classification with the MGMT-pyro  
240 kit.

## 241 **Discussion**

242 To the best of our knowledge, this is the first study to examine all 98 sites with the MGMT pro-  
243 moter CpG island, along with it's shores in multiple patient biopsies.

244 We can conclude that nanopore sequencing of the MGMT promoter region performs as well  
245 or better than standard methods such as pyrosequencing. This is true for both cas9 target-  
246 ted sequencing of the MGMT promoter and inclusion of the MGMT promoter into an adaptive  
247 sequencing panel. Distinct subgroups within both methylated and unmethylated samples are  
248 captured via nanopore sequencing, it will be very interesting to see if there is a difference in  
249 patient outcome between these clusters.

Table 1: Summary of samples included in this study.

	<b>DenStem</b>	<b>Radium</b>	<b>Rapid-CNS</b>	<b>Total</b>
<b>Astrocytoma</b>	3	1	3	7
<b>Astrocytoma HG</b>	0	4	4	8
<b>Pilocytic astrocytoma</b>	0	0	4	4
<b>Glioblastoma</b>	13	29	49	91
<b>Meningioma</b>	0	12	0	12
<b>Metastasis</b>	0	7	0	7
<b>Oligodendroglioma</b>	0	2	6	8
<b>Other</b>	0	10	1	11
<b>Total</b>	16	65	67	148

Table 2: Summary of reported optimal cut-offs for determining methylated versus unmethylated samples

Author	Year	Method	Patients	CpGs	Optimal cut-off	Comment	Reference
Hegi	2019	qMSP	4041		>1.27	"Grey-zone" patients benefit from TMZ	[Hegi et al., 2019]
Johannessen	2018	qMSP, PSQ	48		7 %	PSQ gives better results than other methods	[Johannessen et al., 2018]
Nguyen	2021	PSQ	109		21 %	Higher methylation correlates with longer OS	[Nguyen et al., 2021]
Quillien	2012	MSP, PSQ, MS-HRM	100	5	8 %	PSQ performs best	[Quillien et al., 2012]
Xie	2015	PSQ	43		10 %	Not testing cut-off	[Xie et al., 2015]
Yuan	2017	PSQ	84	4	12.50 %	Higher methylation correlates with longer OS	[Yuan et al., 2017]
Brigliadori	2016	PSQ	105	10	30 %	"Grey-zone" patients do not benefit from TMZ	[Brigliadori et al., 2016]
Radke	2019	PSQ, sqMSP	111		10 %	Best results when PSQ and MSP were combined	[Radke et al., 2019]
Chai	2021	PSQ	173	4	10 %	MGMT promoter methylation has predictive value in IDH-mutant glioblastoma	[Choi et al., 2021]
Dovek	2019	qMSP	165		>1	"Grey-zone" patients benefit from TMZ, higher methylation does not correlate with longer OS	[Dovek et al., 2019]
Siller	2021	MSP, Sseq	215	25		Linear correlation between number of methylated CpG sites and OS	[Siller et al., 2021]

Table 3: Patients used in survival analysis

Sample ID	Age	Sex	Diagnosis	IDH	Resection	Treatment	OS (months)	Status	Pyro_state	NP cluster
1701-2275	66	F	GBM	Neg	GTR	Stupp	14.99	Dead	UnMethylated	1
1701-2430	78	M	GBM	Neg	GTR	Stupp	5.19	Dead	Methylated	2
1701-2590	58	M	GBM	Neg	STR	Stupp	24.5	Dead	Methylated	2
1701-2623	57	F	GBM	Neg	STR	Stupp	28.77	Dead	Methylated	2
1701-2769	73	M	GBM	Neg	STR	Stupp	20.91	Dead	UnMethylated	1
1701-2950	77	M	GBM	Neg	STR	Stupp	11.97	Dead	UnMethylated	1
1501-1486	60	M	GBM	Neg	GTR	Stupp	29.26	Dead	Methylated	2
1501-1757	65	M	GBM	Neg	STR	Stupp	29.69	Dead	Methylated	2
1501-1858	62	F	GBM	Neg	STR	Stupp	6.9	Dead	UnMethylated	1
1501-1880	64	M	GBM	Neg	STR	Stupp	25.48	Dead	Methylated	2
1501-2159	58	M	GBM	Neg	STR	Stupp	21.6	Dead	Methylated	2
1501-2348	58	M	GBM	Neg	STR	Stupp	11.44	Dead	UnMethylated	1
1501-2391	72	F	GBM	Neg	STR	Stupp	21.21	Dead	Methylated	2
1501-2425	58	F	GBM	Neg	STR	Stupp	13.61	Dead	Methylated	2
1601-0227	66	M	GBM	Neg	STR	Stupp	21.96	Dead	Methylated	2
1601-0353	51	M	GBM	Neg	GTR	Stupp	12.85	Dead	Methylated	2
T20-061	64	F	GBM	Neg	STR	Stupp	8.3	Dead	Methylated	1
T20-192	52	F	GBM	Neg	STR	Stupp	23	Dead	Methylated	2
T21-173	66	F	GBM	Neg	STR	Stupp	13.6	Dead	Methylated	2
T21-214	49	M	GBM	Neg	GTR	Stupp	9.4	Dead	UnMethylated	1
T21-216	46	F	GBM	Neg	GTR	Stupp	14.31	Alive	Methylated	2
T21-224	60	F	GBM	Neg	GTR	Stupp	14.08	Alive	Methylated	1
T21-240	55	M	GBM	Neg	GTR	Stupp	15.16	Alive	UnMethylated	2
T21-242	66	M	GBM	Neg	GTR	Stupp	14.47	Alive	Methylated	2
T21-326	39	F	GBM	Neg	GTR	Stupp	13.78	Alive	Methylated	1






Figure 1: Organization of the MGMT promoter. MSP refers to the typical primer sites of methylation specific PCR to determine MGMT promoter methylation. Asterixes represent the 4 CpGs analysed by the Qiagen® MGMT pyrosequencing kit.



(b)

Figure 2: Overview of samples and sequence depth. (a) Classification of all samples used in this study, separated by known methylation status (b) Methylated versus unmethylated samples by method of acquisition (Adaptive sampling, multiplexed nCats, single sample nCats). No bias in sequence depth was observed between methylated and unmethylated samples but single sample runs generally have higher sequence depth than barcoded samples or adaptive sampling.






Figure 3: Comparison of nanopore sequencing and Qiagen® Pyrosequencing kit of CpGs 76-79 in exon 1 of the MGMT promoter. Results show per-site methylation percentage of each CpG (a) or average values of the 4 CpG sites analysed by the Qiagen® MGMT Pyro kit. Black horizontal and vertical lines mark the 10 % cut-off value between methylated and unmethylated samples, as determined by pyrosequencing. Comparison of pyrosequencing classification into methylated versus unmethylated based on a 10% average methylation threshold of CpGs 76-79 in the MGMT promoter (c). The Y-axis represents average methylation percentage of the same four CpG sites based on nanopore sequencing.



(a)



(b)

Figure 4: Something about classifying tumor by an algorithm that only uses 2 CpGs (a). Something about how the different datasets classify differently (b). Something about classifying and sub-classifying tumors according to the methylation of the last 25 CpGs in the MGMT promoter region, as was proposed by Siller *et al.* [Siller *et al.*, 2021] (c). The Y-axis represents aggregated methylation of CpGs 74 to 98 by Nanopore sequencing.






Figure 5: Clustered heatmap of all samples based on nanopore sequencing of CpG island of the MGMT promoter. n = 128




Figure 6: Dotplot showing average methylation percentage of CpG sites in and around the MGMT promoter. Grey areas show 95% confidence intervals of regression lines. n = 128




Figure 7: Heatmap showing unsupervised clustering of glioblastoma samples based on nanopore sequencing of the CpG island in the MGMT promoter. n = 78






Figure 8: K-means clustering of glioblastoma samples.



(a) Pyro classification



(b) Nano Classification

Figure 9: Patient survival based on Pyrosequencing classification (a) or Nanopore Sequencing classification (b)

## References

- [Bady et al., 2016] Bady, P., Delorenzi, M., and Hegi, M. E. (2016). Sensitivity Analysis of the MGMT-STP27 Model and Impact of Genetic and Epigenetic Context to Predict the MGMT Methylation Status in Gliomas and Other Tumors. *Journal of Molecular Diagnostics*, 18(3).
- [Bady et al., 2012] Bady, P., Sciuscio, D., Diserens, A. C., Bloch, J., Van Den Bent, M. J., Marosi, C., Dietrich, P. Y., Weller, M., Mariani, L., Heppner, F. L., McDonald, D. R., Lacombe, D., Stupp, R., Delorenzi, M., and Hegi, M. E. (2012). MGMT methylation analysis of glioblastoma on the Infinium methylation BeadChip identifies two distinct CpG regions associated with gene silencing and outcome, yielding a prediction model for comparisons across datasets, tumor grades, and CIMP-status. *Acta neuropathologica*, 124(4):547–560.
- [Brandner et al., 2021] Brandner, S., McAleenan, A., Kelly, C., Spiga, F., Cheng, H. Y., Dawson, S., Schmidt, L., Faulkner, C. L., Wragg, C., Jefferies, S., Higgins, J. P., and Kurian, K. M. (2021). MGMT promoter methylation testing to predict overall survival in people with glioblastoma treated with temozolomide: A comprehensive meta-analysis based on a Cochrane Systematic Review. *Neuro-Oncology*, 23(9):1457–1469.
- [Brigliadori et al., 2016] Brigliadori, G., Foca, F., Dall’Agata, M., Rengucci, C., Melegari, E., Cerasoli, S., Amadori, D., Calistri, D., and Faedi, M. (2016). Defining the cutoff value of MGMT gene promoter methylation and its predictive capacity in glioblastoma. *Journal of Neuro-Oncology*, 128(2).
- [Choi et al., 2021] Choi, H. J., Choi, S. H., You, S. H., Yoo, R. E., Kang, K. M., Yun, T. J., Kim, J. H., Sohn, C. H., Park, C. K., and Park, S. H. (2021). MGMT promoter methylation status in initial and recurrent glioblastoma: Correlation study with DWI and DSC PWI features. *American Journal of Neuroradiology*, 42(5).
- [Christmann et al., 2011] Christmann, M., Verbeek, B., Roos, W. P., and Kaina, B. (2011). O6-Methylguanine-DNA methyltransferase (MGMT) in normal tissues and tumors: Enzyme activity, promoter methylation and immunohistochemistry. *Biochimica et Biophysica Acta (BBA) - Reviews on Cancer*, 1816(2):179–190.
- [Dovek et al., 2019] Dovek, L., Nguyen, N. T., Ozer, B. H., Li, N., Elashoff, R. M., Green, R. M., Liao, L., Leia Nghiemphu, P., Cloughesy, T. F., and Lai, A. (2019). Correlation of commercially avail-

able quantitative MGMT (O-6-methylguanine-DNA methyltransferase) promoter methylation scores and GBM patient survival. *Neuro-Oncology Practice*, 6(3).

[Hegi et al., 2019] Hegi, M. E., Genbrugge, E., Gorlia, T., Stupp, R., Gilbert, M. R., Chinot, O. L., Burt Nabors, L., Jones, G., Van Criekinge, W., Straub, J., and Weller, M. (2019). MGMT promoter methylation cutoff with safety margin for selecting glioblastoma patients into trials omitting temozolomide: A pooled analysis of four clinical trials. *Clinical Cancer Research*, 25(6):1809–1816.

[Jain et al., 2016] Jain, M., Olsen, H. E., Paten, B., and Akeson, M. (2016). The Oxford Nanopore MinION: Delivery of nanopore sequencing to the genomics community. *Genome Biology*, 17(1).

[Johannessen et al., 2018] Johannessen, L. E., Brandal, P., Myklebust, T. Ø., Heim, S., Micci, F., and Panagopoulos, I. (2018). MGMT gene promoter methylation status – Assessment of two pyrosequencing kits and three methylation-specific PCR methods for their predictive capacity in glioblastomas. *Cancer Genomics and Proteomics*, 15(6):437–446.

[Laver et al., 2015] Laver, T., Harrison, J., O'Neill, P. A., Moore, K., Farbos, A., Paszkiewicz, K., and Studholme, D. J. (2015). Assessing the performance of the Oxford Nanopore Technologies MinION. *Biomolecular Detection and Quantification*, 3.

[Nguyen et al., 2021] Nguyen, N., Redfield, J., Ballo, M., Michael, M., Sorenson, J., Dibaba, D., Wan, J., Ramos, G. D., and Pandey, M. (2021). Identifying the optimal cutoff point for MGMT promoter methylation status in glioblastoma. *CNS Oncology*, 10(3).

[Ostrom et al., 2020] Ostrom, Q. T., Patil, N., Cioffi, G., Waite, K., Kruchko, C., and Barnholtz-Sloan, J. S. (2020). CBTRUS statistical report: Primary brain and other central nervous system tumors diagnosed in the United States in 2013–2017. *Neuro-Oncology*, 22(Supplement\_1):IV1–IV96.

[Patel et al., 2022] Patel, A., Dogan, H., Payne, A., Krause, E., Sievers, P., Schoebe, N., Schrimpf, D., Blume, C., Stichel, D., Holmes, N., Euskirchen, P., Hench, J., Frank, S., Rosenstiel-Goidts, V., Ratliff, M., Etminan, N., Unterberg, A., Dieterich, C., Herold-Mende, C., Pfister, S. M., Wick, W., Loose, M., von Deimling, A., Sill, M., Jones, D. T., Schlesner, M., and Sahm, F. (2022). Rapid-CNS2: rapid comprehensive adaptive nanopore-sequencing of CNS tumors, a proof-of-concept study. *Acta neuropathologica*, 143(5):609–612.

[Quillien et al., 2012] Quillien, V., Lavenu, A., Karayan-Tapon, L., Carpentier, C., Labussi re, M., Lesimple, T., Chinot, O., Wager, M., Honnorat, J., Saikali, S., Fina, F., Sanson, M., and Figarella-Branger, D. (2012). Comparative assessment of 5 methods (methylation-specific polymerase chain reaction, methylight, pyrosequencing, methylation-sensitive high-resolution melting, and immunohistochemistry) to analyze O6-methylguanine-DNA- methyltransferase in a series of 100 glioblastoma patients. *Cancer*, 118(17).

[Radke et al., 2019] Radke, J., Koch, A., Pritsch, F., Schumann, E., Misch, M., Hempt, C., Lenz, K., L bel, F., Paschereit, F., Heppner, F. L., Vajkoczy, P., Koll, R., and Onken, J. (2019). Predictive MGMT status in a homogeneous cohort of IDH wildtype glioblastoma patients. *Acta neuropathologica communications*, 7(1).

[Siller et al., 2021] Siller, S., Lauseker, M., Karschnia, P., Niyazi, M., Eigenbrod, S., Giese, A., and Tonn, J. C. (2021). The number of methylated CpG sites within the MGMT promoter region linearly correlates with outcome in glioblastoma receiving alkylating agents. *Acta neuropathologica communications*, 9(1):35.

[Stupp et al., 2009] Stupp, R., Hegi, M. E., Mason, W. P., van den Bent, M. J., Taphoorn, M. J., Janzer, R. C., Ludwin, S. K., Allgeier, A., Fisher, B., Belanger, K., Hau, P., Brandes, A. A., Gijtenbeek, J., Marosi, C., Vecht, C. J., Mokhtari, K., Wesseling, P., Villa, S., Eisenhauer, E., Gorlia, T., Weller, M., Lacombe, D., Cairncross, J. G., and Mirimanoff, R. O. (2009). Effects of radiotherapy with concomitant and adjuvant temozolomide versus radiotherapy alone on survival in glioblastoma in a randomised phase III study: 5-year analysis of the EORTC-NCIC trial. *The Lancet Oncology*, 10(5).

[Stupp et al., 2017] Stupp, R., Taillibert, S., Kanner, A., Read, W., Steinberg, D. M., Lhermitte, B., Toms, S., Idbaih, A., Ahluwalia, M. S., Fink, K., Di Meco, F., Lieberman, F., Zhu, J. J., Stragliotto, G., Tran, D. D., Brem, S., Hottinger, A. F., Kirson, E. D., Lavy-Shahaf, G., Weinberg, U., Kim, C. Y., Paek, S. H., Nicholas, G., Burna, J., Hirte, H., Weller, M., Palti, Y., Hegi, M. E., and Ram, Z. (2017). Effect of tumor-treating fields plus maintenance temozolomide vs maintenance temozolomide alone on survival in patients with glioblastoma a randomized clinical trial. *JAMA - Journal of the American Medical Association*, 318(23).

[Wongsurawat et al., 2020] Wongsurawat, T., Jenjaroenpun, P., De Loose, A., Alkam, D., Ussery, D. W., Nookaew, I., Leung, Y. K., Ho, S. M., Day, J. D., and Rodriguez, A. (2020). A novel

341 Cas9-targeted long-read assay for simultaneous detection of IDH1/2 mutations and clinically  
 342 relevant MGMT methylation in fresh biopsies of diffuse glioma. *Acta Neuropathologica*  
 343 *Communications*, 8(1):1-13.

344 [Xie et al., 2015] Xie, H., Tubbs, R., and Yang, B. (2015). Detection of MGMT promoter methyla-  
 345 tion in glioblastoma using pyrosequencing. *International Journal of Clinical and Experimental*  
 346 *Pathology*, 8(2).

347 [Yuan et al., 2017] Yuan, G., Niu, L., Zhang, Y., Wang, X., Ma, K., Yin, H., Dai, J., Zhou, W., and Pan,  
 348 Y. (2017). Defining optimal cutoff value of MGMT promoter methylation by ROC analysis for  
 349 clinical setting in glioblastoma patients. *Journal of Neuro-Oncology*, 133(1).

350 [Zhang et al., 2011] Zhang, J., F.G. Stevens, M., and D. Bradshaw, T. (2011). Temozolomide:  
 351 Mechanisms of Action, Repair and Resistance. *Current Molecular Pharmacology*, 5(1).